

AQUADRY

DIGITAL **FITNESS** SCALE

- ✓ *Glass Platform Scale*
- ✓ *Kg / Lb / St Units*
- ✓ *Auto Power Off*

AQUADRY

UP TO
180Kg

INSTRUCTION MANUAL

RELIABLE READING in seconds. Step on the auto-calibrated platform and obtain accurate measurements in pounds / kilograms / stones; weigh up to 400 lb (180 kg).

PRECISION BUILT with 4 high-precision sensors, delivering results with 0.2 lb / 0.1 kg accuracy, clearly displayed on an easy-to-read LCD screen.

TEMPERED STRENGTH Built with a thick and durable 6 mm tempered glass platform enough to provide years of fitness tracking.

EASY TO USE AND STORE Immediate reading as soon as you step on it. Its slim and sleek design saves space. Track your progress with convenience.

ANTI-SKID PADDING Using 4 anti-slip and scratch-resistant pads to keep the weighing platform stable while taking measurements, prevent slipping and protect your floor.

Specifications:

Capacity:	400lb / 180kg
Accuracy:	0.2 lb / 0.1kg
Units:	lb / kg /st
Display:	LCD: 75 x 30mm
Surface Material:	6mm Tempered Glass
Dimensions:	30 x 30 x 2 cm
Item Weight:	2 kg
Battery Powered:	1 x 3V CR2032 Lithium Battery

FEATURES:

- Step-on Technology
- Auto On/Off
- Auto Calibration
- Lb/Kg/St Unit Conversion
- Slim and Sleek Design
- Anti-Skid Padding
- Low Battery Indicator: Lo
- Overloading Indicator: Err

HOW TO OPERATE:

1. Pull the plastic protection from the battery to initialize the scale.
2. Set the preferred unit, Lb/Kg/St at the bottom part of the scale.
3. For best accuracy, place the scale on a hard, flat surface. Avoid carpet or soft surfaces.
4. Step on the scale gently to turn on and read your weight. Relax and stand still until stable reading is displayed.
5. It will automatically turn off when not in use.

NOTES:

1. "Lo" symbol is displayed when low battery
2. "Err" symbol is displayed when overloaded

CARE AND MAINTENANCE:

1. Remove the battery if the scale will not be used for a long time.
2. Clean the scale with a soft cloth moistened with water or a mild detergent. Never use gas, hot water or other harsh chemicals to clean the scale.
3. Avoid extreme temperature, humidity, direct sunlight, moisture and dust.
4. Never shock, beat, drop or attempt to open the scale and tamper the mechanism or circuitry.
5. The scale is only designed for personal home use. It's not for professional use in hospitals, medical facilities and trade.
6. The included battery is intended for demonstration purposes only. We recommend replacing a new battery.

CORRECT DISPOSAL OF THIS PRODUCT



This marking indicates that this product should not be disposed of with other household wastes. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmentally safe recycling.

AQUADRY

DIGITAL
FITNESS
SCALE

MADE UNDER LICENSE IN CHINA

Sourced & distributed exclusively for the Philippines by:
ASIA PACIFIC COSMETICS PHILS. CORPORATION
9/F SMB Bldg., No. 40 San Buenaventura St., Bagoing Ilog,
Pasig City, Philippines 1600

Warranty:

This product is unconditionally guaranteed for one year
against all defects of workmanship and materials. This guarantee
does not affect your statutory rights.

DESIGNED IN



FOR FURTHER ASSISTANCE OR INFORMATION

CALL OUR CUSTOMER SERVICE ON:

Telephone No: (+632) 8465 9000 or 8465 9094 Mobile No: (+63) 923 746 5611

Email Add: customerservice@primeline.com.ph

 www.aquadry.com.ph  [@aquadryph](https://www.facebook.com/aquadryph)  [aquadry_ph](https://www.instagram.com/aquadry_ph)